Symbols Used on Labeling
Auf der Verpackung verwendete Symbole
Symboles utilisés sur l'étiquetage
Simboli usati sull'etichetta
Símbolos utilizados en las etiquetas
Símbolos utilizados na embalagem

REF

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English

A. REFERENCES

- These instructions were developed using the guidance given in the following standard:
- ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
 - AAMI TIR12: Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: Richtlinie für Produkthersteller.

B. CLEANING AND STERILIZATION PROCEDURE

Assembly/disassembly of these devices is required. When properly performed, cleaning and/or sterilization do not compromise the use and mechanical performance of these devices. These devices are used with or on patients who may harbor both recognized and unrecognized infections. To prevent the spread of infection, all reusable devices must be thoroughly cleaned and sterilized after use on each patient.

C. LIMITATIONS ON REPROCESSING

Repeated processing does affect the useful life of these devices. End of useful life (which is typically 1 year from the start of use) is normally determined by wear due to repeated processing, wear due to use and / or damage.

D. VALIDATION

The following recommended sterilization methods have been validated in compliance with federal and international guidance standards. Cleaning and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility/end user to perform the appropriate validation testing.

E. CONTAINMENT AND TRANSPORTATION

It is recommended that devices are reprocessed as soon as is reasonably practical following use.

F. CAUTIONS

- Never immerse the device in water or cleaning detergent.
- Never clean the device with bleach, chlorine-based detergents, liquid or chemical disinfectants, or any products containing sodium hydroxide (such as, INSTRUKLEZ or Buell Cleaner). These fluids will degrade the anodized aluminum coating and potentially damage the handpiece's functionality.
- Never clean the device in an ultrasonic cleaner or a combination washer/sanitizer.
- To avoid damaging the handpiece, do not impact or subject the handpiece to blunt force. Rough handling of the handpiece may cause damage to the handpiece.
- Do not use the Arthrex handpiece for any purpose other than its intended use. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- Only use the APS II console that has been developed by Arthrex specifically for the handpieces.
- Do not attach the handpiece during APS II Console Self Test or the Brightness or Language Selection Programming Modes. Refer to the *APS II User's Guide* for more information.
- To avoid damaging the shaver handpiece and disposables used with the shaver handpiece, do not subject the disposable to excessive force. Excessive force may render the disposable and/or shaver handpiece inoperable.

G. PREPARATION FOR CLEANING

- Ensure that the soaker cap is securely covering the connector end of the cable.
- Open all suction pathways to enable maximum exposure. For devices equipped with a removable suction valve, remove it by depressing the center of the pivot face and sliding out the valve.
- Remove dried-on soil from devices, especially in areas such as joints and crevices, prior to washing.

H. MANUAL CLEANING

- Perform rinsing and cleaning as soon as possible after use.
- Thoroughly scrub the device with a clean, soft brush that is dampened with a mild, enzymatic, pH-balanced detergent. Always avoid any harsh materials that can scratch or mar the surface of the device.
- Remove all traces of blood, debris, and stains.
- Do not immerse the device in soap solution or rinse-water.
- To clean the cannulated section of the device:
 - Feed the wire end of a cleaning brush through the cannulation of the device.
 - Occulate the brush back and forth from both the front and the rear of the device. Brush all sections of the cannulation.
 - It is recommended that the cannulation be inspected with an illuminated, magnifying scope to ensure that all debris is removed.

Note: After the cannulated section of the device and the suction valve (when removable) are clean, reinsert the suction valve by sliding it into position until an audible click is heard.

- Manipulate all moving parts of the device to ensure that all soil is removed. If any soil remains, clean the device again until the soil is removed.
- Keep the nose of the device pointed downward and rinse it under running water to remove all traces of soap.
- Flush the cannulation thoroughly with water, keeping the nose pointed down during the flushing operation.

I. INSPECTION AND MAINTENANCE

- Arthrex shaver handpieces are precision medical instruments and must be used with care.
- Inspect the device for damage prior to use and at all stages of handling thereafter.
- If damage is detected, consult your Arthrex representative for guidance.
- Check device for visible soil. It is recommended that the cannulation be inspected with an illuminated, magnifying scope. Clean the device using the guidelines for manual cleaning if any soil is visible.
- Lubricate all moving parts with a soluble instrument lubricant prior to sterilization.**

J. STERILIZATION

- Sterilization does not eliminate the need for proper cleaning of devices. Sterilization capabilities, cleaning, handling, and storage of instrumentation are the responsibility of qualified facility/user personnel.
- The APS II Handcontrol Shaver Handpiece may be sterilized by using one of the following autoclaving parameters:

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	4 Minutes	20 to 30 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	132°C – 135°C (270°F – 275°F)	10 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

CAUTION: NEVER allow the APS II console plug pins to have any contact with liquids. Remove dust or moisture, if necessary, with dry compressed air. ONLY dry connectors may be plugged into the APS II console.
Cooling – The device must be adequately cooled, after being removed from the sterilizer. It should not be touched during the cooling process. Do not place the device on a cold surface or immerse in a cold fluid.

DFU-0154r10

CAUTION: NEVER submerge the cable or handpiece in liquid to cool, as this can damage the cable, handpiece, connectors, and/or seals.

K. PACKAGING

- A standard packaging material may be used. Ensure that the pack is large enough to contain the device.
- Handpieces may be loaded into dedicated instrument trays or general-purpose sterilization trays. Ensure that the handpiece is protected and does not exceed 8.5 kg/18.7 lb per tray. Wrap the trays using the appropriate method.

L. STORAGE

Store in a cool, dry location.

M. SPECIAL PRECAUTION - TRANSMISSIBLE SPONGIFORM ENCEPHALOPATHY AGENTS

- It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy (TSE) Agents.
- The agents for transmission of Creutzfeldt-Jakob disease are believed to be resistant to normal processes of disinfection and sterilization and therefore the normal processing methods of decontamination and sterilization as described above may not be appropriate where CJD transmission is a risk.
- In general, the issues that come into contact with orthopedic surgical instruments are those of low TSE infectivity. However, particular precautions should be taken when handling instruments that have been used on known, suspected, or at-risk patients.

N. DIRECTIONS FOR USE

Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

1. Foot Control Option

- Insert the connector from the handpiece into the handpiece receptacle of the console so that the red dots on the connector and receptacle align and engage easily.
- The shaver handpiece cable connects and locks to the console to prevent accidental separation during use. To avoid damage, only disconnect the shaver handpiece by pulling on the cable connector.
- The console will recognize the handpiece type, and display available options for the operator to adjust, as desired.
- The sterile suction hose is pushed onto the suction connection.
- Adjust the flow volume with the suction control.
 - Move the suction control switch to the "+" on the handpiece body (fully opened) to obtain the maximum flow volume.
 - Move the suction control switch to the "-" on the handpiece body (fully closed) to stop flow.
 - The suction can be adjusted between the above settings to achieve appropriate level.
- To insert a shaver blade or burr, fully insert the hub into the hub bore of the handpiece. Watch the locking pin slide back and forth to check that it is locked in place. There will be an audible "click" sound when the handpiece is securely locked.
- The locking pin of the shaver handpiece is compatible with all multi-use and disposable Arthrex blades and burs.
- To remove a shaver blade or burr, slide the exposed locking pin to the opposite side and remove the blade.
- Use the RPM selection button on front panel of the console, or the selection mechanism on the handpiece, the PushPad, or the footswitch to adjust the speed of the handpiece.
- The shaver handpiece supports speeds from 500 rpm to 8,000 rpm, depending on the direction of operation. Refer to the *APS II User's Guide* for more details.
- The Reverse, Forward, and Oscillate buttons on the APS II console select which direction the shaver handpiece motor rotates. Refer to the *APS II User's Guide* for more information.

II. Hand Control Option – AR-8325H, AR-8330H and AR-8330RH

NOTE: APS II console must have software revision 4.31 or greater for the Handcontrol handpiece to operate.

- To activate the AR-8325H handpiece, press and release the button. To deactivate the handpiece, press and release the button a second time. To activate the AR-8330H or AR-8330RH handpiece, press the button. To deactivate the handpiece, press the button a second time.
- To change direction, press and release the button.
- To lower the speed of the AR-8325H handpiece while it is in operation, press and hold the button. An audible tone will be heard from the console, indicating that the speed is being altered.
- To increase the speed of the AR-8325H handpiece, press and hold the button. An audible tone will be heard from the console, indicating that the speed is being altered.
- Note:** While in the speed change mode, you can alternate between the and the buttons to adjust the speed accordingly. Normal button operation will resume when neither button is depressed for 3 seconds and another audible tone is heard from the console.
- To change the speed of the AR-8330H or AR-8330RH handpiece, press and hold the button until an audible tone is heard from the console, indicating that the speed may be altered. Continue pressing the button to increase speed. Press the button to decrease speed. Normal button operation will resume when neither button is pressed for 3 seconds and another audible tone is heard from the console.
- To transfer control of the buttons on the Handcontrol handpiece to the footswitch, either connect the footswitch:
 - to the same channel as that of the Handcontrol handpiece, or
 - to the other channel and then press Toggle on the footswitch.

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