

NEW 10" INSPEKTOR: 3X Larger Image




INSPEKTOR® CT Video Borescope Inspection System




The Waterproof Video Borescope for use in GI and SPD Lumen Inspection


- A Plug-N-Play System with 10" Monitor - Computer Not Required
- Automatic White Balance and Illumination – Saves Time
- May be High Level Disinfected and is STERRAD Compatible
- Small Enough to get into Bronchoscopes and Long Enough for GI Endoscopes



Dirty Arthroscopic Shaver



Debris Left in Endoscope



Dirty Duodenoscope

Two Sizes

INSPEKTOR CT ●
2.1 mm x 1.8 m

INSPEKTOR CT MINI ●
1.2 mm x 1.0 m