


# Multi-Center Flexible Lumen Inspection Survey

**70%**  
of Inspected Instruments Were Dirty

Participating  
Hospitals : **71**

Instruments Inspected with a 1mm or 2mm  
Waterproof Lumen Inspection System : **860**

It is hard to clean  
what you cannot  
see. You can guess  
it is clean, but will  
not know until it is  
inspected



To be sterile,  
instruments have  
to be as clean on  
the inside as well  
as the outside.



## Instrument Types With The Highest Percentage of Foreign Matter in the Lumens

Arthroscopic Shavers Hand Pieces	81% (76 out of 94)
Suction Tubes	73% (348 out of 480)
Cannulated Orthopedic Instruments	72% (101 out of 141)
Laposcopic Instruments	67% (75 out of 112)

